

Blue Whiz - Urinal and Bowl Cleaner

PRODUCT #: PK-4710

CHARACTERISTICS:

Appearance: Blue Liquid
 Odor: Characteristic
 pH: N/A
 Specific Gravity: 1.048
 Weight: 8.74
 Viscosity: Water Like
 Flash Point: N/A
 Dilution: Ready To Use



ADVANTAGES:

Blue Whiz is a fast-acting, effective, easy to use professional-grade bowl and urinal cleaner. It makes a whiz of your bowl cleaning biz!

DIRECTIONS:

FOR USE IN URINALS AND AS A TOILET BOWL CLEANER: Thoroughly saturate bowl mop with cleanser. Rub bowl until discolorations and lime deposits disappear. If badly discolored, repeat above and allow to stand for 5 minutes before flushing.
FOR CLEANING THE INVISIBLE TRAP OR OUTLET IN BOWLS: Pour fluid into bowl, leave for several hours or overnight. Insoluble incrustations that have accumulated in the trap will be removed. To assure clean traps and bowls, repeat at least once a week.

PRECAUTIONS:

Inhalation: Remove affected person to fresh air. If breathing difficulty persists, seek medical attention
Skin: Remove contaminated clothing; wash affected area with soap and water; if irritation persists, seek medical attention. Launder contaminated clothing before reuse.
Eyes: Check for and remove any contact lenses. Immediately flush eyes with clear running water for 15 minutes while holding eyelids open. Seek medical attention immediately.
Ingestion: Drink two glasses of water followed by milk, milk of magnesia or other non-alcoholic liquids; DO NOT induce vomiting; seek medical attention immediately.

TRANSPORT INFO:

Proper Shipping Name - Corrosive liquids, N.O.S.,
 (Contains Hydrochloric Acid)
 DOT Hazard Class - 8
 Packing Group - III
 ID# - UN 1760



AVAILABLE IN:

Gal. Jug, 5 Gal. Bucket, 55 Gal. Drum,
 275 Gal. Tote, 300 Gal. Leg Tank or delivered via Tanker Truck



GUARANTEE: Jamson Labs' / Power Kleen's manufacturing and laboratory controls ensure uniform quality throughout our product line. All Jamson Labs / Power Kleen manufactured products are guaranteed for 1 year when used as directed.